
SORIN | SYNTHESIS 4D PH.I.S.I.O.

GB – ENGLISH - DIRECTION FOR USE.....	5
IT – ITALIANO – ISTRUZIONI PER L’USO	12
FR - FRANÇAIS – MODE D’EMPLOI	20
DE – DEUTSCH – GEBRAUCHSANWEISUNG INHALT	28
ES – ESPAÑOL – INSTRUCCIONES DE USO	36

011654/00 – 03/2009


Fig. 1

Reference Fig. 1

1. Venous Reservoir
2. Venous Inlet
3. Venous Temperature Probe Fitting
4. Quick prime port
5. Reservoir Outlet (to Centrifugal Pump)
6. Venous/arterial sampling stopcock manifold
7. Vent/vacuum port
8. Filtered cardiomy inlets
9. Filtered cardiomy Luer ports
10. Over-under pressure safety valve
11. Reservoir mount locking ring
12. Centrifugal pump inlet (from reservoir)
13. Arterial outlet
14. Gas Inlet
15. Gas Scavenging connector
16. Blood Cardioplegia outlet
17. Arterial temperature probe fitting
18. Water inlet and outlet connectors
19. Arterial filter
20. Recirculation/purge stopcock
21. Recirculation/purge line
22. Centrifugal pump
23. Arterial Pressure Monitoring Port
24. Arterial sampling port

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Fig. 2 - Soporte Synthesis 4D

TABLE OF CONTENTS

A. DESCRIPTION
B. TECHNICAL FEATURES
C. INTENDED USE
D. SAFETY INFORMATION
E. PREPARATION AND SETUP
F. PRIMING AND RECIRCULATION PROCEDURE
G. STARTING PERFUSION
H. DURING PERFUSION
I. TERMINATING PERFUSION
J. BLOOD RECOVERY AFTER PERFUSION
K. DEVICE CHANGE-OUT
L. USE OF ACTIVE VENOUS DRAINAGE WITH VACUUM
M. POSTOPERATIVE DRAINAGE
N. MEDICAL DEVICES FOR USE WITH SYNTHESIS 4D
O. RETURN OF USED PRODUCTS
P. LIMITED WARRANTY

A. DESCRIPTION

SYNTHESIS 4D is an integrated system for hemodynamic support, gas exchange in extracorporeal perfusion and cardiopulmonary bypass for periods of up to 6 hours. SYNTHESIS 4D comprises the following elements integrated into one assembly: a microporous hollow-fibre membrane oxygenator and heat exchanger, an arterial filter, a centrifugal pump and a hard-shell venous reservoir.

The arterial filter includes a 40 µm screen. The centrifugal pump is equipped with a vane impeller designed to move blood by centrifugal force.

SYNTHESIS 4D is coated with Phosphorylcholine (Ph.I.S.I.O.). Devices coated with Ph.I.S.I.O. are used when a coated blood path is desired. The Ph.I.S.I.O. coating improves the blood compatibility of the device by reducing platelet adhesion to the coated surfaces.

The device is single-use, non-toxic, non-pyrogenic, supplied **STERILE** and individually packaged, sterilised by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

The device is available in the following versions:

[A] SYNTHESIS 4D with integrated arterial filter, centrifugal pump and venous reservoir for adult bypass.

B. TECHNICAL FEATURES

Maximum flow rate	8 l/min
Membrane type	Microporous polypropylene
Membrane surface area	2.0 m ²
Heat exchanger surface area	0.14 m ²
Arterial filter:	
- Filter screen	40 µm
- Type	polyester
- Surface area	0.04 m ²
Static priming volume (post use recovered)	451 ml
Priming volume	525 ml
(oxygenator module + heat exchanger + arterial filter + centrifugal pump)	
Venous Reservoir:	
- Max. Cardiotomy flow rate	4 l/min
- Maximum reservoir volume	4300 ml
- Minimum reservoir volume	300 ml
- O/U pressure relief valve	open at P<-80 mmHg (-10.4kPa/-0.10 bar/-1.53psi) and at P>5 mmHg (0.7 kPa/ 0.007 bar/ 0.1 psi)
Connections :	
- Venous inlet	1/2" (12.7 mm)
- Arterial filter outlet	3/8" (9.53 mm)
- Blood Cardioplegia outlet	pos-lock
Maximum Pressure Rating:	
- Blood Pathway	750 mmHg (100 kPa/ 1bar / 14,5psi)
- Gas Pathway at outlet	0 mmHg (0 kPa/ 0 bar/ 0 psi)
- Water Pathway	2250 mmHg (300 kPa/ 3 bar / 43,5 psi)

C. INTENDED USE

SYNTHESIS 4D is intended to be used in adult surgical procedures requiring extracorporeal gas exchange support, to control the arterial/venous temperature, as a venous blood reservoir /filter element to eliminate gas emboli and remove blood component aggregates larger than 40 µm.

SYNTHESIS 4D must not be used longer than 6 hours. Contact with blood for longer periods is inadvisable.

SYNTHESIS 4D is intended for use only with the Sorin Group Deutschland Stöckert Centrifugal Pump Console.

SYNTHESIS 4D must be used in combination with the medical devices listed in paragraph N (Medical devices for use with SYNTHESIS 4D).

D. SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

⚠ WARNING

WARNING indicates serious adverse reactions and potential safety hazards for the practitioner and/or the patient that may occur in proper use or misuse of the device and also the limitations of use and the measures to be adopted in such cases.

⚠ CAUTION

CAUTION indicates any special care to be exercised by a practitioner for safe and effective use of the device.

EXPLANATION OF SYMBOLS USED ON THE LABELS

 For single use only (Do not reuse)

 Batch code (number)
(reference for product traceability)

 Use by (Expiry date)

 Date of manufacture

 Sterile - Ethylene oxide sterilised

 Non Pyrogenic

 Contains PHTHALATE

 Latex free

 Warning: Do not resterilize.

 Contents sterile only if package is not opened, damaged or broken

 Catalogue (code) number

 Attention, see instruction for use

 This way up

 Fragile; handle with care

 Keep away from heat

 Keep dry

ea Units

The following is general safety information with the aim of advising the operator in preparing to use the device.

⚠ WARNING

- The user should carefully check the device for leaks during set-up and priming. Do not use if any leak is detected.
- The device must be used in accordance with the instructions for use provided in this manual.
- Also specific information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

- For use by professionally trained personnel only.
- SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.
- Do not use the device if it has been cracked, dropped or otherwise physically damaged.
- FRAGILE. Handle with care.
- Keep dry. Store at room temperature.
- Always apply and maintain a correct dose and accurate monitoring of the anticoagulant before, during and after the bypass.
- Do not use the SYNTHESIS 4D for a period longer than 6 hours. The system has not been qualified through in vivo, ex vivo, or clinical studies for long term use (longer than 6 hours) as a bridge to transplant, or for pending recovery of the natural heart or for extracorporeal circulation (ECMO). The use of the device may result in pump failure, reduced priming capability, excessive blood trauma, degradation and/or corrosion of blood contacting materials with the consequent possibility of particles passing through the bypass circuit to the patient, leaks, and increased potential for gaseous emboli entering the arterial line.
- For single use and for single patient use only: during use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body and due to its specific design it cannot be fully cleaned and disinfected at the end of use. Therefore, reuse on other patients might cause cross-contamination, infection and sepsis. In addition, reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device must not undergo any further processing.
- Do not resterilize.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if STERILE.
- The device contains phthalates. Considering the nature of body contact, the limited contact duration and the number of treatments per each patient, the amount of phthalates which might be released from the device do not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.
- The water temperature at the heat exchanger inlet must not exceed 42°C (108 °F).
- The water pressure in the heat exchanger must not exceed 2250 mmHg (300 kPa / 3 bar / 44 psi).
- Do not allow the blood side pressure to exceed 750 mmHg (1 bar/ 100 kPa / 14.5 psi).
- Do not allow the pressure on the blood side of the membrane to fall below the pressure on the gas side of the membrane. Blood compartment pressure must always be maintained higher than gas compartment pressure to prevent formation of gas emboli in the blood compartment.
- SYNTHESIS 4D should always be placed at a level not higher than the patient.
- Ensure that the flow of the main pump is always greater than that of the cardioplegia or haemoconcentration pump. In any event, the flow of the cardioplegia line must never exceed 1 l/min.
- Avoid occluding, even accidentally, of the gas outlet line: any over-pressurization of the gas compartment may generate micro-emboli in the blood compartment.
- Check the patency of the gas circuit by circulating gas in the circuit. If the gas circuit is occluded, replace the oxygenator.
- Arrange the surgical cannula in such a way that air is not aspirated into the venous line.
- Do not ventilate the SYNTHESIS 4D while the arterial and venous lines are clamped, as this can cause evaporation of prime solution across the microporous membrane and result in air in the arterial line.
- Prevent that halogenated liquids such as halothane and fluothane come into contact with the polycarbonate structure of the device. This would cause damage which would compromise the integrity and functionality of the device.
- When using the SYNTHESIS 4D venous reservoir as postoperative thoracic drainage, carefully follow the instructions in paragraph M "POSTOPERATIVE DRAINAGE" provided in this document.
- Do not for any reason separate the gas module/heat exchanger/arterial filter/pump from the venous reservoir. Only perform this operation if the venous reservoir is to be used for postoperative thoracic drainage and only when extracorporeal circulation has been terminated.
- Do not separate the SYNTHESIS 4D modules in order to place the gas module at a level or position different from its original configuration at the time of opening the sterile package.
- When using the SYNTHESIS 4D venous reservoir applying the

method of active venous drainage with vacuum, carefully follow the instructions in paragraph L "USE OF ACTIVE VENOUS DRAINAGE WITH VACUUM" in this user manual.

- The special positive and negative overpressure valve mounted on the SYNTHESIS 4D venous reservoir assures perfect postoperative functioning of the system: it starts operating and releases positive pressures above 5 mmHg (0.7 kPa / 0.007 bar/0.1 psi) and negative pressures below -80 mmHg (-10,4 kPa / -0.10 bar/-1,53 psi). Do not for any reason occlude the external access hole of pressure relief valve. Remove tab inserted in pressure relief valve from the valve prior to use.
- The vent cap must be removed in order to prevent inadvertent pressurization of the reservoir.
- The vent must be open at all times during the operation of the reservoir or be attached to a regulated vacuum source not to exceed -50 mmHg (-6,66 kPa / -0,07 bar / -0,97 psi) in procedures using vacuum assisted venous return.
- Obstruction of the vent port may result in pressurization of the reservoir and potentially gaseous air bubbles passing air emboli to the patient and/or damage to the device.
- Do not allow the venous reservoir to become pressurized as this could obstruct venous drainage, force air retrograde into the patient, or cause air to enter the blood path of the oxygenator.
- SYNTHESIS 4D may only be used with the appropriate HOLDER.
- For further information and/or in case of complaints contact SORIN GROUP ITALIA or the authorised local representative.

CAUTION

- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
- Do not operate the pump unprimed. Doing so may damage the pump.
- To aid in the prevention of gas or particulate embolization, SORIN GROUP ITALIA recommends the use of safety devices, including bubble sensors and pre-bypass filters, for all procedures using the SYNTHESIS 4D.
- Entry of massive air into the pump will cause the device to deprime and the blood flow to stop.
- The pump must not be run unattended.
- Inner surfaces of the system are Ph.I.S.I.O. coated, currently SORIN GROUP ITALIA is not aware of any contraindication to use of systems having components treated with Ph.I.S.I.O.

E. PREPARATION AND SETUP

1) POSITIONING THE HOLDER

Position the SYNTHESIS 4D HOLDER on a pole and fix it by means of the clamp behind the load bearing "arch." Be sure to orient the holder such that the device will be upright, and check that there is sufficient space for the pump driver.

2) FIXING THE DEVICE ON THE HOLDER

WARNING

- Do not use if sterile packaging is damaged, unsealed, or has been exposed to moisture or other conditions that would compromise the sterility of the device.
- Check the expiry date on the label. Do not use the device after this date.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.

Remove the device from the sterile packaging.

WARNING

- Visually inspect and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Particularly, turn the device upside down and verify the integrity of the four fins connecting the arterial filter to the oxygenating module.
- Do not use solvents such as alcohol, ether, acetone, etc., as contact may cause damage to the device.
- Do not allow halogenated liquids such as Halothane and Fluothane to come into contact with the polycarbonate housing of the device. This could cause damage which may compromise the integrity and proper functioning of the device.

Fix the SYNTHESIS 4D oxygenator onto the holder. Check that the water connector locking lever is in the "OPEN" position. Make sure that the water connectors of the oxygenator are aligned with the Hansen connectors of the holder. Insert the oxygenator Hansen connectors and push SYNTHESIS 4D down into the holder.

Turn the locking lever to the "CLOSED" position.

Slide the pump driver into place and secure it to the pump head of SYNTHESIS 4D. Verify that the centrifugal pump is correctly engaged with the pump driver locking tab. Close the support collar around the

pump driver.

Only when the locking lever shows "CLOSED" the oxygenator is correctly inserted in the holder.

⚠ WARNING

- **Manually verify that the arterial filter and oxygenating module are stably connected.**

3) PRIMING THE THERMOCIRCULATOR

Connect the water tubes to the holder by means of the female Hansen connectors (SORIN GROUP ITALIA code 09028 or equivalent).

⚠ WARNING

- **Connect the water tubes coming from the thermocirculator to the oxygenator using Sorin Group Italia female Hansen connectors.**
- **The use of different connectors from those indicated may cause resistance inside the circuit and reduce the efficiency of the heat exchanger.**
- **Do not obstruct the hole in the lower cover on the heat exchanger as it is the outlet of the safety channel which helps prevent fluids from crossing from one compartment to another.**

4) CHECKING THE HEAT EXCHANGER

The device is checked by recirculating water in the exchanger for a few minutes. The integrity of the structure is guaranteed if there are absolutely no water leaks from the compartment or from the hole of the safety channel.

5) DEVICE SET-UP AND CONNECTIONS

⚠ WARNING

- **All connections downstream of the pump must be secured by means of ties.**

Once the SYNTHESIS 4D has been positioned in the holder and the centrifugal pump driver is put in place, complete the setup of the remaining extracorporeal circuit components.

- VENOUS RETURN. A 1/2" (12.7 mm) venous line can be attached to the connector marked on the venous reservoir as "VENOUS RETURN LINE" (Ref.2, Fig.1). The venous return connector can be rotated 360° as needed.
- SUCTION LINES. After removing the protective caps from the "filtered" inlets on the top of the venous reservoir (four 1/4" / 6.35 mm inlets, two 3/8" / 9.53 mm inlets), connect the ends of the suction tubes (Ref. 8, Fig.1) and rotate the turret orienting the filtered inlets towards the suction pumps as desired.
- ARTERIAL LINE: Verify that the 3/8" (9.5 mm) arterial line is connected to the arterial filter outlet port. Connect the centrifugal pump flow probe to circuit where desired.(Ref. 13 , Fig.1)
- Remove the yellow cap from the vent connector marked "VENT/VACUUM PORT" (Ref. 7, Fig.1) of the venous reservoir.
- Verify that the arterial recirculation/purge line (Ref. 21, Fig.1) is connected to the recirculation/purge Luer of the arterial filter, and connect the opposite end of the recirculation/purge line to the reservoir.

⚠ CAUTION

- **If oxygenated blood for blood cardioplegia needs to be used, remove the red pos lock and connect the 1/4" blood line of the cardioplegia circuit to the arterial blood access port of SYNTHESIS 4D (Ref. 16, Fig.1) by means of the adaptor (provided with the product).**

⚠ WARNING

- **If using the blood cardioplegia outlet during bypass, the line to be connected must be open and not under pressure so that the blood can flow into it.**

6) SAMPLING SYSTEM

- If it is preferred to take samples from a more distant position, remove the A/V sampling manifold (fig. 1, ref. 6) from its housing and attach it on the dedicated holder. The tubing lines of the sampling system allow positioning the manifold in a range of about 1 meter. Position the stopcock handles toward the access ports of the manifold (A/V shunt position).
- Check the connection of the arterial sampling line to the Luer next to the recirculation/purge Luer at the top of the arterial filter.

⚠ WARNING

- **The arterial sampling site connector Luer does not include a one-way valve. Verify that there is a one-way valve in the sampling line to prevent accidental air introduction into the arterial line.**
- Check the connection of the venous sampling line to a Luer port provided on the venous side of the device.

⚠ CAUTION

- **Check that all the Luer connections are secure. All the**

accessory lines connected to the device must be connected tightly in order to prevent accidental introduction of air into the device or blood loss.

7) TEMPERATURE AND PRESSURE MONITORING

Use SORIN GROUP ITALIA temperature probes (code 09026) or equivalent.

- Connect the arterial temperature probe to the temperature probe fitting (Ref. 17, Fig. 1), and connect the venous temperature probe to the temperature probe fitting provided on the venous side of the device (Ref. 3, Fig. 1).
- To measure the oxygenator inlet pressure, connect a pressure measurement line to the venous Luer close to the heat exchanger inlet.
- To measure oxygenator outlet pressure, connect a pressure measurement line to the Luer port (Ref. 23, Fig. 1) close to the arterial blood access port.

8) VENTILATING GAS SUPPLY

Remove the green cap from the connector marked "GAS INLET" (Fig. 1, Ref. 14) and connect the 1/4" gas line. The gas must be supplied through a suitable air/oxygen mixer. The connection for a capnograph is obtained on the central axis of the "GAS ESCAPE" connector (Fig. 1, Ref. 15).

⚠ CAUTION

- **The "gas escape" system was designed to prevent any potential risk of occlusion of the gas outlet; such occlusion would cause immediate passage of air to the blood compartment.**

9) VAPOROUS ANAESTHETICS

The oxygenator is suitable for single use with volatile anaesthetic isoflurane and sevoflurane, by means of a suitable narcosis gas evaporator.

If these vaporous anaesthetics are used, some method of scavenging the gas from the oxygenator should be considered. The protocol, the concentration and monitoring of the anaesthetic gases administered to the patient is under the sole responsibility of the physician performing the treatment.

An integral 3/8" (9.5 mm) gas scavenging fitting is available for this purpose. It is located in the centre of the gas escape system. Connect a piece of 3/8 x 3/32" (9.5 x 2.4 mm) tubing to this fitting and to a negative pressure source with a flow rate in excess of the ventilating gas flow rate.

⚠ WARNING

- **The only volatile anaesthetics suitable for this use are isoflurane and sevoflurane.**
- **The methods adopted for vaporous anaesthetic gas scavenging should not in any way increase or reduce the pressure level at the oxygenator fibres.**

F. PRIMING AND RECIRCULATION PROCEDURE

⚠ WARNING

- **Do not use alcohol based priming solutions as they can compromise the structural integrity and functionality of the device.**

1) CHECK THE HEAT EXCHANGER

Circulate water through the heat exchanger prior to priming the oxygenator. Check the water system for leaks, proper water temperature range, and adequate water flow. If a water to blood leak exists, water will collect in the tubing connected to the oxygenator inlet/centrifugal pump outlet port or drip from the hole in the bottom of the oxygenator module indicating the device must not be used. The maximum pressure rating of the heat exchanger is 2275 mmHg (300 kPa/ 44psi).

2) KEEP THE GAS FLOW CLOSED

3) CHECK THAT THE RECIRCULATION/PURGE LINE IS OPEN

4) OCCLUDING THE VENOUS AND ARTERIAL LINES

5) VENOUS RESERVOIR PRIMING

Secure with ties all the aspiration lines connected to the venous reservoir. Fill the venous reservoir with sufficient liquid to ensure that the intended hematocrit is obtained, taking into account that:

- the static priming volume of the oxygenator already including the arterial filter is 525 ml;
- the 3/8" tube contains 72 ml/m;
- the 1/2" tube contains 127 ml/m.

6) PRIMING THE CIRCUIT

⚠ WARNING

- **The pressure in the blood compartment of the oxygenator**

module must not exceed 750 mmHg (100 kPa / 1 bar / 14 psi).

- In order to maintain positive hydrostatic pressure during the priming and recirculation phase, keep the reservoir above the oxygenator.
- It is recommended to use a pre-bypass filter to capture any particles that may be present in the circuit or in the priming solution.
- a. While not required, flushing the circuit with carbon dioxide gas will aid in priming.
- b. Open the venous line attached to the outlet of the reservoir and slowly allow priming solution into the SYNTHESIS 4D.
- c. Allow the priming solution to prime the pump, heat exchanger, oxygenator and arterial filter.
- d. Verify that the connection tube between the venous reservoir outlet and the inlet connector of the centrifugal pump is completely primed and debubbled.
- e. Activate centrifugal pump to recirculate priming solution through recirculation/purge line of the arterial filter.

WARNING

- During this phase with the recirculation/purge line open and the arterial line clamped, never exceed a pump rate of 1000 rpm.
- f. Remove the clamp from the arterial line, increase the flow rate and recirculate the priming fluid through the A/V loop in order to prime and debubble the circuit.
- g. Set the pump flow rate to 6 litres/min. Continue to recirculate for 3-5 minutes. While recirculating, check the system for air bubbles and tap the entire circuit to facilitate the removal of air.
- h. Ensure that the sampling manifold system is primed. The sample system will prime itself during recirculation when the stopcocks are positioned to allow flow through the sample lines and manifold. The handles point to the off position. Flow through the sample system should be stopped when there is not flow in the bypass circuit.
- i. Ensure that all other parts of the circuit are primed and de-aired.
- j. Once the circuit and device are primed, stop the centrifugal pump and clamp the arterial and venous lines. Close the recirculation/purge stopcock.
- k. Prior to bypass, disconnect the A/V loop, clamp the arterial and venous lines, then pass the clamped lines to the surgeon for connection to the respective cannulae.

CAUTION

- It is intended that systemic anticoagulation be utilized while the pump is in use. Anticoagulation levels should be determined and verified by the physician based on risks and benefits to the patient and monitored throughout the case.

WARNING

- If the cardioplegia circuit has been connected to the arterial blood access port, check that the circuit has been primed.
- If tubing lines used for pressure measurement are connected, check that the lines have been primed.
- Occlude the line with a clamp some centimetres away from the outlet.
- Do not apply negative pressures to the coronary outlet. Negative pressures in the blood compartment may cause the formation of gas micro-emboli.

G. STARTING PERFUSION

1) INITIATING FLOW

Start the centrifugal pump.

WARNING

- Do not operate the pump for extended periods of time if the arterial outlet to the patient is clamped and the purge/recirculation stopcock is closed at the same time. The temperature inside the pump may increase, increasing the risk of cellular damage.
- Do not operate the pump with the centrifugal pump inlet clamped. A negative pressure would be generated in the pump and air bubbles could form.
- Do not remove the arterial outlet clamp at the beginning of extracorporeal circulation until an adequate outlet pressure has been reached in order to prevent backflow. The pump must be operated to provide a pressure higher than the patient systemic pressure and the circuit head pressure. Monitor the pump flow rate, its rpm and the systemic pressure as indices of potential backflow.
- Carefully monitor for signs of occlusion throughout the circuit.
- Do not remove the arterial outlet clamp at the beginning of extracorporeal circulation until an adequate outlet pressure has been reached in order to prevent backflow. The pump must be operated to provide a pressure higher than the patient systemic pressure and the circuit head pressure. Monitor the pump flow rate, its rpm and the systemic pressure as indices of potential backflow.

2) INITIATING BYPASS

WARNING

- Control the flow rate by adjusting the pump speed. Attempting to adjust the blood flow rate by partial clamping of the device outlet may augment the risk of cellular damage.
- When stopping circulation, clamp the device outlet and stop pump rotation immediately, or reduce the pump speed to less than 1000 rpm if the arterial filter recirculation line is open.

First remove the clamp from the arterial line and then from the venous line. Enter into bypass with a blood flow suited to the size of the patient. Constantly check the blood level in the venous reservoir.

3) ESTABLISH GAS FLOW

The suggested gas/blood flow ratio in normothermia is 1:1 with an FiO₂ of 80:100%.

WARNING

- Always activate the gas flow after the blood flow. The gas/blood flow ratio must never exceed 2:1.

4) FUNCTIONING OF HEAT EXCHANGER

Check the temperature of the venous and arterial blood.

5) MANAGING THE ARTERIAL FILTER

SYNTHESIS 4D has an integrated 40 µm self-debubbling filter.

The filter is equipped with a stopcock with two connections, one for the re-circulation/purge line, and the other for the arterial sampling line. With the stopcock of the re-circulation/purge line in the automatic purging position, the air can continuously be evacuated from the oxygenator module during bypass.

It is recommended to continuously evacuate the air for the first 10-15 min. of extracorporeal circulation.

6) BLOOD GAS ANALYSIS

After a few minutes of bypass, the amount of gas in the blood must be checked. In consideration of the values read, act as follows:

High pO ₂	→	decrease FiO ₂
Low pO ₂	→	Increase FiO ₂
High pCO ₂	→	Increase gas flow
Low pCO ₂	→	Decrease gas flow

Close the re-circulation/purge stopcock.

H. DURING PERFUSION

1) CHECKING THE VENOUS RETURN

If a higher venous return is required, lower the oxygenator/venous reservoir level with respect to the patient.

WARNING

- The ACT (Activated Coagulation Time) must always be greater than or equal to 480 seconds in order to assure proper anticoagulation of the extracorporeal circuit.
- If it is necessary to administer anticoagulant to the patient, use the A/V sampling stopcock.
- The minimum operating volume of the cardiectomy reservoir is 300 ml. However, to ensure an adequate response time in case of venous inflow obstruction, it is recommended to maintain an adequate volume in addition to the 300 ml minimum level. Do not exceed 4300 ml in the venous reservoir.
- Avoid that in the heat exchanger the difference between the water temperature and the blood temperature exceeds 15°C. Higher values may cause the formation of bubbles from the gases dissolved in the blood

2) ARTERIAL SAMPLING

Insert a sample syringe into the arterial sample stopcock Luer. Position the stopcock handles of the manifold to allow arterial blood to flow through the manifold. The pressure on the arterial side will allow flow. Draw the sample of blood from the arterial sample stopcock. Turn off the arterial stopcock before removing the syringe.

WARNING

- The arterial sampling site is under pressure. Always insert a syringe into the sampling site before opening the stopcock.

3) VENOUS SAMPLING

Ensure that the arterial stopcock is closed. Insert a sample syringe into the venous stopcock Luer. Open the central stopcock and draw at least 10-15 ml of blood prior to taking the venous sample. Close the central and the venous stopcock. Return this blood through one of the filtered Luer connections positioned on the top of the reservoir. Open the venous stopcock and draw a sample of venous blood and close the stopcock before removing the syringe.

4) DRUG INJECTION

Insert the medication syringe into the Luer connector of the central stopcock. Open the central and venous stopcocks and inject the drug

into the manifold and venous sample line.

Close the central stopcock to the flow of the medication syringe and allow an arterial-venous "wash" through the stopcock manifold. Turn the stopcocks to the closed position when "washing" has been completed.

CAUTION

- **Take the blood from the stopcocks only when the pump is running. If not, the pressure in the blood compartment would decrease causing the formation of air bubbles.**

5) LOW-FLOW RECIRCULATION OFF PUMP

(hypothermia associated with circulatory arrest).

- Reduce the gas flow to a value below 500 ml/min. and reduce FiO_2 to a value not exceeding 30% in order to prevent reaching high non-physiological pO_2 values.
- Open the re-circulation/purge line and occlude the venous reservoir inlet line.
- Occlude the arterial line of the oxygenator module.
- Re-circulate at 1000 ml/min for the entire duration of the circulatory arrest.
- To re-enter into bypass from the circulatory arrest, open the venous and arterial lines and slowly increase the blood flow.
- Set the re-circulation/purge line to the automatic purging position.
- Align the gas flow and FiO_2 depending on the needs of the patient.

WARNING

During the venous reservoir emptying phase at particularly low levels and/or reduced flows, greater care must be taken.

CAUTION

If during a particularly long ECC, condensation forms in the gas outlet with a corresponding drop of pO_2 and an increase of pCO_2 , the performance of the oxygenator could be improved by a short but substantial increase in the gas flow. For example, "washing" of 10 sec. with a gas flow of 2 l/min could be sufficient. Should this not prove effective, do not repeat washing, but replace the oxygenator.

I. TERMINATING PERFUSION

It should be terminated depending on the conditions of each individual patient. Act as follows:

- Stop flow through the sample system prior to concluding bypass.
- Conclude bypass according to normal operating procedures.
- Turn the gas flow off.
- Slowly reduce the centrifugal pump speed while gradually occluding the arterial line.
- Open the recirculation/purge line to the recirculation position.
- Occlude the venous line between the patient and the reservoir.
- Increase the pump speed to 1000 ml/min maximum.
- Recirculate as required.

CAUTION

- **If extracorporeal circulation has to be restored, maintain a minimum blood flow inside the oxygenator (maximum 1000 rpm).**
- **Do not turn the thermocirculator off during the re-circulation phase.**
- **Check that the cardioplegia circuit connected to the coronary outlet is properly occluded.**

J. BLOOD RECOVERY AFTER PERFUSION

The volume of blood inside the SYNTHESIS 4D can be recovered according to normal operating procedures. Use aseptic techniques in handling the device during volume recovery.

K. DEVICE CHANGE-OUT

A spare oxygenator/pump must always be available during perfusion. After 6 hours of use with blood or if situations arise such that, on judgment of the person responsible for perfusion, the safety conditions for the patient are compromised (insufficient oxygenator performance, leaks, anomalous blood parameters etc.), proceed as follows to replace SYNTHESIS 4D.

Supplies needed for the change-out procedure:

- One (1) sterile tubing cutter
- Eight (8) tubing clamps
- One (1) sterile SYNTHESIS 4D
- One (1) SYNTHESIS 4D holder

CAUTION

- **Use a sterile technique throughout the replacement procedure.**

- 1) Connect a second oxygenator/pump assembly holder to the pump console. Place the holder as close to the first holder as is possible, such that the new oxygenator/pump will be oriented in the same orientation as the old oxygenator.
- 2) Remove the new device from packaging and inspect sterile wrapper for damage. Do not use if sterile barrier or device has been compromised or damaged. Place the new device in the second holder and mount the pump driver per § E2. The inlet and outlet ports of the new device should be oriented in the same direction as the ports of the device being replaced. Stop the flow through the sample system of the old device, and verify that the recirculation/purge valve is closed and that the arterial access line is clamped.
- 3) Connect the water lines to the new device. Circulate water through the heat exchanger and check for leaks in the tubing connected to the blood inlet port to the hole on the bottom of the oxygenator.
- 4) Terminate bypass according to normal procedures. This should involve clamping of the venous line and arterial line.
- 5) Turn off ventilating gas. Remove the gas line from the old oxygenator module and place it on the new oxygenator module.
- 6) Connect the arterial sample line of the new sample system kit to the arterial sampling port on the arterial filter. Remove the old venous sample line from the venous sampling port of the venous reservoir, and attach the venous sample line of the new sample system. Make sure the stopcock handles on the sample manifold are positioned to allow flow through the sample lines and manifold while ensuring there is no leakage from the sample site ports.
- 7) Connect the recirculation/purge line of the new device to the reservoir.
- 8) Place two tubing clamps on the 3/8" (9.5 mm) pump inlet line. This is the line coming from the reservoir to the inlet port of the old device. Place these clamps close to the pump inlet port and approximately 1.5 inches (3.8 cm) apart from each other. Using aseptic technique, cut the line between the clamps. Cut near the clamp closest to the pump.
- 9) Without removing the clamp, attach the 3/8" (9.5 mm) pump inlet line (coming from the reservoir) to the new pump inlet port.
- 10) Place two tubing clamps on the 3/8" (9.5 mm) filter outlet line. Place these clamps close to the outlet port and approximately 1.5 inches (3.8 cm) apart from each other. Using aseptic technique, cut the line between the clamps. Cut near the clamp closest to the filter outlet port.
- 11) Without removing the clamp, attach the 3/8" (9.5 mm) outlet line to the new filter outlet port.
- 12) Disconnect the arterial blood access line from the old device and reconnect to the arterial access port on the new device.
- 13) Remove the clamp on the pump inlet line of the new device.
- 14) Open the recirculation/purge valve to the purge position.
- 15) Check the security of all the connections. Ensure that there is still a clamp on the arterial line to the patient.
- 16) Prime the new device and evacuate any air using the priming and recirculating procedure using the AV bypass loop in the circuit.
- 17) Turn off the sample system. Turn off the pump and close the recirculation/purge line.
- 18) Remove the clamps on the venous and arterial lines, and reinstate bypass according to normal procedures.
- 19) Remove the arterial temperature probe from the old device and connect to the arterial temperature probe fitting of the new device.
- 20) Remove the gas scavenging line from the old device and connect to the new device.

L. USE OF ACTIVE VENOUS DRAINAGE WITH VACUUM

This method may be applied at any time during extra-corporeal circulation, provided that the use of the device follows the constraints described below. Using an active venous drainage tubing kit supplied separately and a vacuum regulation device, SYNTHESIS 4D can be used in active venous drainage applications with vacuum assist. Set up the regulation device according to the manufacturer's instructions.

- 1) Aseptically open the packaging of the kit for active venous drainage with vacuum.
- 2) Connect the end with the blue cap to the vent connector of the venous reservoir, marked "VENT/VACUUM PORT", and the end with the red cap to a vacuum regulating device. The latter must be connected to a vacuum source.
- 3) Close the clamp and the green cap on the line connected to the reservoir.
- 4) If considered necessary to interrupt venous assist, remove the yellow cap and open the clamp on the line described in 3).

WARNING

- **Use aseptic technique and sterile tubing for the above mentioned steps. Secure tubing ports by sliding tubing over each port until fully engaged. Use of a locking tie is**

recommended on all tubing connections. Ensure the cap on any unused port is airtight prior to initiating vacuum assisted venous return.

- It is advisable not to exceed -50 mmHg (-6.66kPa / -0.07 bar / -0.97 psi) negative pressure applied to the reservoir.
- Periodically check functioning of the vacuum regulating device and the degree of vacuum.
- Utilization of vacuum assisted venous return can lead to negative pressures in the oxygenator and the potential for air to be pulled across the oxygenator membrane into the blood pathway. The sample system, the arterial purge line, a haemoconcentrator, and the centrifugal pump may also provide a conduit for the vacuum to be applied to the arterial side of the oxygenator.
- The use of the following is recommended in procedures utilizing vacuum assisted venous return: reliable wall suction and calibrated vacuum regulator, vapour trap in vacuum line, venous line pressure monitor, vacuum relief valve in vacuum line, one-way valves where necessary to prevent retrograde flow.

M. POSTOPERATIVE DRAINAGE

The SYNTHESIS 4D reservoir is also intended to be used for postoperative chest drainage collection and transfusion.

CONTRAINDICATIONS

The SYNTHESIS 4D reservoir is contraindicated for autotransfusion use under the following circumstances:

- When pericardial mediastinal, pulmonary or system infection exists;
- When malignancy, gross contamination, lymphatic failure, or a perforated intestine exists or is feared to exist;
- After a return to surgery for exploration (not contraindicated if massive bleeding occurs); and
- When protamine is used before the reservoir is removed from the extracorporeal circuit.

WARNING

- Before use ensure that the SYNTHESIS 4D venous reservoir has not been polluted as a result of handling which has compromised sterility.
- The sterility of the system depends on the degree of asepsis resulting from use of the venous reservoir and the care taken in connecting the kit during the conversion to thoracic drainage.
- Ensure the cap on any unused port is airtight prior to initiating chest drainage collection.
- Complications such as coagulopathies, particulate and air embolism, and blood coagulation may occur during postoperative autotransfusion. The physician is responsible for determining the quality of the blood product before it is returned to the patient.
- For postoperative use, never clamp the chest drainage tubes, unless only momentarily to change out the device. Blockage of air or fluid movement can lead to tension pneumothorax or subcutaneous emphysema.
- Never autotransfuse blood that has resided in the reservoir for more than six (6) hours. This will reduce the potential for bacterial contamination. After autotransfusion has been stopped, the system can continue to be used for chest drainage collection until the chest drainage tubes are removed from the patient.
- If, during postoperative use, blood is ever observed in the filter bypass line, the subsequent autotransfused blood should first be passed through a microfilter to remove any clots that may be present. Replacement of the reservoir in this situation is recommended.
- Maintain a minimum of 50 ml of blood in the reservoir during autotransfusion. If this is not done, gaseous bubbles could be passed to the patient.
- During autotransfusion use an infusion pump with a bubble detector to avoid the possibility of gaseous bubbles being passed to the patient.
- Prepare the vacuum regulating device with integrated valve for use, following the instructions of the manufacturer provided with the device.
- Should negative pressure be applied to the system such that the overpressure valve intervenes, the recovered blood will no longer be suitable for transfusion.
- Periodically check functioning of the vacuum regulating device.
- if the recovered blood exceeds the amount established by postoperative routine, immediately notify the physician.
- If the recovered blood is to be re-infused, an increase in vacuum will be noted in the reservoir. A filtered venting system must therefore be arranged in order to maintain the original negative operating pressure.
- Check the degree of contamination of the blood collected related to the time of use of the drainage system. The results of these

tests must be taken into account prior to blood transfusion.

- The recovered blood must be re-transfused with a suitable volumetric infusion pump or with the syringe and three-way stopcock technique used manually as an alternative to the pump.
- If the reservoir used does not have an integrated microfilter, fit the re-infusion line with a 25-40 micron blood filter.

CONVERTING THE SYNTHESIS 4D RESERVOIR FROM INTRAOPERATIVE TO POSTOPERATIVE USE

Supplies needed for conversion:

- One (1) reservoir holder
 - One (1) sterile water manometer/seal
 - One (1) sterile 1/4" (6.4 mm) water manometer line
 - One (1) sterile 1/2" (12.7 mm) ID occluding plug
 - One (1) sterile IV pump line adapter
 - One (1) sterile 3/8" (9.5 mm) chest drain line
 - One (1) foot (30 cm) of sterile 1/4" (6.4 mm) ID tubing
 - Two (2) sterile 3/8" (9.5 mm) sealing caps
 - Four (4) sterile 1/4" (6.4 mm) sealing caps
1. Ensure that the extracorporeal circuit will not be needed for an emergency return to use for bypass support.
 2. Remove the venous inlet tubing from the SYNTHESIS 4D reservoir. Using aseptic technique, seal the venous inlet port with the 1/2" (12.7 mm) ID occluding plug.
 3. Using aseptic technique, connect the 1/4" (6.4 mm) ID tubing from a filtered to an unfiltered 1/4" (6.4 mm) port such that the tubing loops vertically above the reservoir for maximum visibility. Blood in this tubing will indicate the filter has become occluded and is being bypassed.
 4. One by one, using aseptic technique, carefully remove the tubing from the remaining ports used intraoperatively and plug the ports with sealing caps.
 5. Ensure the SYNTHESIS 4D reservoir is completely empty. Remove the 3/8" (9.5 mm) ID tubing from the reservoir outlet and replace with the IV pump line adapter.
 6. Mount the venous reservoir holder in a suitable location for patient transport.
 7. Remove the venous reservoir from the oxygenator using the following steps:
 8. Press the tab on the Reservoir Mount Locking Ring to unlock the ring and rotate it clockwise until it stops.
 9. Lift the reservoir away from the oxygenator.
 10. Place the SYNTHESIS 4D reservoir in the dedicated reservoir holder. Place the water manometer/seal on the holder. Attach the 1/4" (6.4 mm) water manometer line to the vent/vacuum port of reservoir and the "liner lid tubing" port of the water seal.
 11. Set up the water seal according to the manufacturer's instructions.
 12. Using aseptic technique, connect the chest drain line(s) to the filtered 3/8" (9.5 mm) port(s) on the reservoir cover. If one of these ports is not used, the cap provided with the reservoir will provide an adequate seal.
 13. The SYNTHESIS 4D reservoir is now ready for transport to the ICU.
 14. SETUP OF SYNTHESIS 4D RESERVOIR FOR AUTOTRANSFUSION IN THE ICU
 15. Once the unit has been moved to the ICU, connect the water manometer/seal to a vacuum source according to the water manometer/seal instructions.
 16. After the initial drainage rate has been determined, and when there is at least 75 ml of volume in the reservoir; reinfusion can begin if indicated.
 17. If blood is to be reinfused, prime an IV pump set by first clamping off the vacuum source line. Connect the inlet of the IV pump line to the outlet port of the IV pump line adapter. Prime the set as recommended by the pump manufacturer following directions for cellular fluids. Install the set into the pump. Remove the clamp from the vacuum source. Connect the distal end of the IV set into a patient IV access site used for volume administration.
 18. Monitor and record the patient's chest drainage rate per your institution's protocol. To minimize the potential for air embolism, a minimum of 50 ml of blood should remain in the reservoir after reinfusion.
 19. If rapid bleeding occurs, the blood may clot and occlude the reservoir's filter. If the filter is occluded, blood may back up and pass through bypass loop and flow directly into the unfiltered section of the reservoir, bypassing the filter. In this situation, use of a blood transfusion filter in the IV line is required during reinfusion.

WARNING

- Use of the SYNTHESIS venous reservoir as thoracic drainage is the responsibility of the intensive care staff.
- Its preparation is the joint responsibility of the person performing perfusion and the intensive care staff.

- Its connection to the patient is the responsibility of the attending surgeon.
- Use by untrained and unqualified personnel is not permitted.
- All the connections to the Sorin reservoir should be carried out with the utmost care and rapidity in order to prevent contamination.
- The connections and the caps should be checked to assure that the system is tight.
- Always keep the device in a vertical position: do not incline it, not even when the patient is being moved.
- The SYNTHESIS reservoir must always be placed at a subthoracic level so that its drainage functions can properly be performed.
- Activate the line vacuum following the instructions provided by the manufacturer of the vacuum regulator.
- Any re-infusion of recovered blood is on the responsibility and at the discretion of the attending physician

N. MEDICAL DEVICES FOR USE WITH SYNTHESIS 4D

The centrifugal pump is intended for use only with Sorin Group Deutschland Stöckert Centrifugal Pump Consoles. Refer to the console operator's manual for console operating procedures.

- SYNTHESIS 4D holder code 050148
 - Reservoir holder code 050546
 - Stopcock holder code 050134
 - Active Drainage and Vacuum Assist tubing packs (as required)
- The circuit connections must be made with tubes of a diameter compatible with the dimensions of the connectors located on the device (3/8", 1/4", 1/2").

Temperature control must be carried out with SORIN GROUP ITALIA probes, code 09026 or YSI Series 400 compatible.

Use a Sechrist air/oxygen blender (Sorin Group Italia code 09046) or a system with compatible technical features to control the blood gas concentration.

Any heating/cooling system (thermocirculator) may be used, provided that the connectors to the oxygenator are the Hansen type (Sorin Group Italia code 09028).

O. RETURN OF USED PRODUCT

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

Sorin Group Italia reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

CAUTION

It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood-borne infectious diseases.

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned goods authorisation (RGA) number is required from Sorin Group USA Inc. prior to shipping.

If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the USA is:

Sorin Group USA, Inc.

Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031.

P. LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended.

SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging.

However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected.

SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device.

SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser.

The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary.

The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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SORIN GROUP USA, Inc